




HOW TO USE: Print this template on A4 paper, landscape. DO NOT SCALE. Print at 100% or disable page scaling. Cut along the dashed lines and remove inner shape. Insert your cat's nose and mouth into the front and side cut-outs to see which mask will fit best.

AeroKat 

Feline Aerosol Chamber



The exclusive **Flow-Vu*** Inhalation Indicator moves toward the cat as they inhale and will only move if there is a good seal.



Mask Sizer

Each **AeroKat*** Chamber is shipped with a small and medium mask. The mask selected should cover the cat's nose and mouth but not the eyes.

Helping animals around the world breathe easier



Trudell Medical International*

www.aerokat.co.uk

BreathEazy[®]
respiratory relief for cats, dogs and horses

154 Worcester Road, Malvern,
Worcestershire WR14 1AA, United Kingdom
T: +44 (0) 845 680 8975 F: +44 (0) 845 680 8976
www.breatheasy.co.uk